

Wireless connectivity accessories

Everything you need  
to **stay connected**



**oticon**  
life-changing technology

# Make everyday tasks easier

Oticon wireless accessories are a comprehensive and easy-to-use range of devices developed to improve your listening and communication experiences with your hearing aids.

Our wireless connectivity accessories help you get the most out of your Oticon Bluetooth® hearing aids. For example, you can stream the sound from your TV, smartphone, music player, PC or other devices to your hearing aids - just like when you wear wireless headphones.

## Control your hearing aids with Oticon ON

The Oticon ON app lets you discreetly adjust the volume, switch between programmes, check battery level, use HearingFitness™, access the Internet of Things, and much more.

The app also lets you control your streaming from multiple TV Adapters, and now offers a streaming equalizer for a personalised listening experience.

## Connect to the world

EduMic helps when you need direct and clear sound, by transmitting audio or a speaker's voice into your hearing aids. The combination of EduMic and Oticon Bluetooth hearing aids gives you a clear signal with exceptional sound quality.





# Stream audio directly from a mobile device

Oticon Bluetooth hearing aids offer an immersive streaming experience with excellent sound quality from mobile devices.

## iPhone®, iPad® and iPod touch®

Oticon Bluetooth hearing aids are Made for iPhone (MFi) hearing aids. They can directly connect to iPhone, iPad and iPod touch for streaming audio and thereby act as wireless stereo headphones – without the need for an intermediary device.

## Android™ devices

Oticon More also supports Audio Streaming for Hearing Aids (ASHA) and can therefore stream audio directly from Android devices that also support ASHA. Users of other Oticon Bluetooth hearing aids or of devices that do not support ASHA need ConnectClip as an intermediary device.



Made for

iPhone | iPad | iPod

Works with

android

Apple, the Apple logo, iPhone, iPad, and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android, Google Play, and the Google Play logo are trademarks of Google LLC.



## Hands-free calls for any modern smartphone

With ConnectClip you can have hands-free calls on an iPhone, Android or any other smartphone. ConnectClip streams the caller's voice directly to your hearing aids while the built-in microphone clearly transmits what you say to the caller. ConnectClip transforms your hearing aids into a high-quality wireless headset enabling you to:


- Make and receive hands-free calls in the car or on the go
- Wirelessly stream music, videos or other sound from your smartphone to your hearing aids
- Stream the voice of your conversation partner to your hearing aids
- Control the volume and change programmes remotely

Works with  
any modern  
smartphone





Hear the person  
talking with  
**more clarity**



Even with the best hearing aids, it can sometimes be difficult to clearly hear what someone is saying – such as in noisy environments, when attending a lecture, talking on your home phone or over the computer.

ConnectClip can help you in all these situations by streaming sound wirelessly and directly to your hearing aids. This will make it easier to listen to the person you want to hear while further reducing the background noise.

### **Hear speech at a distance or in noise**


- Use ConnectClip as a remote microphone by handing it to the person you want to hear better
- ConnectClip streams the speaker's voice directly to your hearing aids
- Use the Oticon ON app to control the speech volume and the level of background noise
- Range: up to 20 metres

### **Chat over the computer**

- Pair ConnectClip directly to your computer\*
- Stream any sound wirelessly from your computer to your hearing aids
- Enjoy video conversations and have your voice streamed back to the computer
- Range: up to 10 metres

### **Talk on your landline phone**

- Connect Phone Adapter 2.0 to your existing phone system
- Use ConnectClip to answer the call
- Enjoy hands-free calls and get clear sound directly in your ears
- Range: up to 20 metres



\*If your computer does not have Bluetooth technology, you need a USB dongle to connect your computer to ConnectClip



## Control your hearing aids easily and discreetly

If you want to control the volume, change programmes or simply mute your hearing aids, you can use any of the three remote control options: ConnectClip, Remote Control or the Oticon ON app.

## Get TV sound delivered straight to your hearing aids

TV Adapter 3.0 lets you enjoy your favourite TV programme or listen to music from your stereo directly in your hearing aids - just as if you were wearing wireless stereo headphones.

### With TV Adapter you can:


- Clearly hear the TV or listen to music by streaming the sound directly to your hearing aids
- Choose the level you want to hear while leaving the TV at a volume that is comfortable for others
- Pair up to four TV Adapters to your hearing aids\*
- Use the Oticon ON app to control the volume of the streamed sound remotely, and switch between multiple TV Adapters
- Range: up to 15 metres



\*All Oticon Bluetooth hearing aids except Oticon Opn and Siya



# Boost your hearing aids with Oticon ON



The Oticon ON app is available for iPhone and Android smartphones, and opens up a wide range of features and functionalities:

- Switch between programmes
- Adjust volume
- Check battery status
- Locate your hearing aids
- Access to Tinnitus SoundSupport™ feature\*
- Use HearingFitness™
- Connect to the Internet of Things
- Control connectivity devices such as ConnectClip, EduMic and multiple TV Adapters
- Fine-tune the sound when streaming audio for a personalised listening experience

## Boost your hearing when you need it most

In an everyday situation where you need more help to make speech stand out from noise, you can use the sound booster function in the Oticon ON app to give your hearing aids an extra boost.

For information on compatibility, please visit:  
**[oticon.global/compatibility](https://oticon.global/compatibility)**

**Read more on [oticon.global/app](https://oticon.global/app)**

\*if enabled by your hearing care professional during the fitting of your hearing aids



## Fine-tune your streaming experience

The new streaming equalizer in the Oticon ON app lets you fine-tune the sound when streaming music or videos. This enables you to create a more personalised listening experience that suits your preferences.

### Getting started

1. Download the free Oticon ON app from the App Store® or Google Play™
2. Pair your Oticon hearing aids with your iPhone, iPad, Android smartphone or tablet
3. You're all set!



Download on the  
**App Store**



GET IT ON  
**Google Play**

Apple, the Apple logo, iPhone, iPad, and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android, Google Play, and the Google Play logo are trademarks of Google LLC.

## Enjoy **remote visits** with your hearing care professional

Although nothing can replace a face-to-face consultation for more major issues, modern technology makes it easy to carry out routine adjustments to your hearing aids without the need for you to go into the clinic.

### The Oticon RemoteCare app enables you to:

- Get expert hearing aid follow-up consultations in the comfort of your own home - or wherever else is convenient\*
- Save time and money on travelling to the hearing care clinic for minor hearing aid adjustments
- See, hear, talk and text with your hearing care professional during your remote visit
- Get your hearing aids adjusted remotely to assist you in getting the best out of them
- Test your new hearing aid settings while you are talking with your hearing care professional



### Download Oticon RemoteCare from the App Store or Google Play

The availability of this service may vary. Please contact your local hearing care professional for more information.

For information on compatibility, please visit: [oticon.global/compatibility](https://oticon.global/compatibility)

\*For information on technical requirements, please visit: [oticon.global/remotecare-support](https://oticon.global/remotecare-support)



## Maximizing the opportunities in life

Hearing key information can be a challenge. Sometimes it's when trying to hear someone speak over background noise or in a crowd. Other times, it's hard to hear at a distance, such as in church, at the theatre, during education, and when playing sport.

EduMic is designed to solve these issues by transmitting an individual voice directly into the user's hearing aids. It is a remote microphone that's easy to use. What's more, it features a telecoil to pick up audio from public audio loop systems and transmits it to the hearing aids.

EduMic also streams audio from all sorts of audio devices like computers, smart boards, tablets or smartphones - connecting via a standard 3.5 mm headphone jack.

### Helping your child hear in noisy situations

Children with hearing loss can be especially challenged in school, where noisy listening environments can make it hard to follow the class. This is where EduMic offers great learning support by giving children direct access to the teacher's voice.

Outside school hours, EduMic can be your invaluable link to your child. At home, it doesn't matter if you have appliances running, the TV on, or siblings bickering - EduMic will bring your voice directly to your child's ears. You can also use EduMic at the playground, out shopping, and in any other noisy environment.



Contact your local hearing care professional:

Oticon wireless connectivity accessories are compatible with the following Oticon hearing aids (all styles except CIC and IIC):

Oticon More™

Oticon Opn S™

Oticon Opn Play™

Oticon Opn™

Oticon Xceed

Oticon Xceed Play

Oticon Ruby

Oticon Siya

For more information on products and compatibility between devices, talk to your hearing care professional today or go to [oticon.global/connectivity](https://www.oticon.global/connectivity)