

Contact Us

6238 8832

www.hearingpartners.com.sg

wecare@hearingpartners.com.sg




SCAN ME



Move Beyond the Noise—Find Relief That Fits You.

Oticon Jet PX is designed for anyone with tinnitus, including those with normal hearing, providing comfort, control, and an improved quality of life.



Ask your hearing care professional about a personalized fitting and discover how easy it is to manage tinnitus on your terms.

*Relief from Tinnitus.
Personalized, discreet,
and always within reach.*



Find Relief from Tinnitus

with a Tinnitus Masker

 Hearing Partners

What is Tinnitus

Tinnitus is humming, ringing, buzzing, or other sounds in the ears or head that lasts more than 5 minutes and occurs at least twice a week.

It's common, even if you have normal hearing, and can disrupt sleep, focus, and daily comfort.

What is a Tinnitus Masker

A tinnitus masker is a device that plays gentle, soothing sounds to help cover up (or "mask") tinnitus.

Unlike commercial sound generators, Oticon Jet PX is a small discreet ear-level device that is specifically designed to manage tinnitus and helping you regain peace of mind.

How Does a Tinnitus Masker Help



By introducing soothing sounds, the masker helps your brain focus less on the tinnitus, making it less noticeable and easier to ignore.



Being an ear-level device, Oticon Jet PX is highly portable and convenient to use—it can be worn at all times for your best comfort.



With regular use, many people find their tinnitus becomes less intrusive, allowing for better relaxation, concentration, and sleep.

What to Take Note Of

- Tinnitus maskers do not cure tinnitus, but they can make it much easier to manage day-to-day.
- Relief is personal—some people notice immediate comfort, while others experience gradual improvement.
- For best results, use your masker consistently and work with a hearing care professional to tailor the sounds to your preference.
- Always use the device at safe volume levels, and follow guidance from your hearing care professional to protect your hearing health.

Oticon Jet PX

Advanced, Discreet, and Designed for You



Slim, stylish, and comfortable for all-day wear, Oticon Jet PX is available in a range of colors to suit your lifestyle.



Powered by Oticon's Deep Neural Network (DNN) and BrainHearing™ technology for natural, balanced sound and optimal listening comfort.



Enjoy direct streaming from your iPhone, iPad, Mac, and select Android devices—stream your favorite sounds or use the tinnitus masker feature for personalized relief.

Available in 5 colours



C092
Steel Grey



C091
Silver Grey



C090
Chroma Beige



C093
Chestnut Brown



C063
Diamond Black

Take Control with Oticon Companion App – ComfortPulse



Personalize Your Relief: The Oticon Companion app lets you adjust your tinnitus masker right from your phone or Apple Watch—no need to touch the device.





ComfortPulse Feature: Choose from constant tones or four unique pulsating sound patterns, and adjust the speed and volume for each ear independently or together.



More Than Masking: Use the app to fine-tune your sound experience, check battery status, and even locate your device if misplaced



Discreet and Convenient: Make changes anytime, anywhere, for instant comfort and peace of mind.



Apple, the Apple logo and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Google Play, and the Google Play logo are trademarks of Google LLC.